



Product Code

| | |
|-----------|---|
| UTGG-1600 | Pyknometer (Bottle Type) 250 ml |
| UTGG-1605 | Pyknometer (Bottle Type) 500 ml |
| UTGG-1610 | Pyknometer (Bottle Type) 1000 ml |
| UTGG-1615 | Pyknometer (Bottle Type) 2000 ml |
| UTGG-1620 | Pyknometer (Bottle Type) 3000 ml |
| UTGG-1625 | Pyknometer (Bottle Type) 5000 ml |
| UTGG-1630 | Double Edged and Capillary Tubed Funnel |

Standards

EN 1097-6

Bottle Type Pyknometers are used to determine the specific gravity of aggregates. 250 ml, 500 ml, 1000 ml, 2000 ml, 3000 ml and 5000 ml.



The Bottle Type Pyknometer is supplied complete with

- Double Edged and Capillary Tubed Funnel

| | Dimensions | Weight (approx.) |
|-----------|-------------|------------------|
| UTGG-1600 | Ø60x330 mm | 0,30 kg |
| UTGG-1605 | Ø80x360 mm | 0,50 kg |
| UTGG-1610 | Ø100x410 mm | 0,60 kg |
| UTGG-1615 | Ø130x450 mm | 0,90 kg |
| UTGG-1620 | Ø160x480 mm | 1 kg |
| UTGG-1625 | Ø180x510 mm | 2 kg |
| UTGG-1630 | Ø40x230 mm | 0,10 kg |